

Product Family



everStick<sup>®</sup>

from GC

fiber reinforcements  
for daily  
dentistry

- Reliable • Easy to use • Minimally invasive • Extra strong
- Aesthetic • Scientifically proven • Cost effective

**everStick<sup>®</sup>**  
GIANT OF FIBERS

**GC**

## everStick®C&B

For minimally invasive fiber reinforced composite bridges

- Surface retained bridges
- Inlay and onlay bridges
- Hybrid bridges and temporary bridges
- Laboratory-made bridges



## everStick®POST

For advanced root canal post and core structures

- Individually formed root canal posts



## everStick®PERIO

For patient friendly splinting

- Periodontal splinting
- Combined periodontal splint and surface-retained bridge



## everStick®NET

For easy and aesthetic splinting of traumatized teeth

- Labial splinting of traumatized teeth
- Labial periodontal splinting
- Repair and reinforcing of veneers




## everStick®ORTHO


For aesthetic retention

- For patient-friendly metal free orthodontic retainer





A cross-sectional view of everStick fiber. Silanated glass fibers are impregnated with PMMA and bis-GMA.



Inside everStick fiber. Individual fibers, bis-GMA and PMMA for the unique IPN structure.

## everStick® for daily dentistry

everStick glass fiber reinforcements have been developed to provide solutions for modern, patient-friendly dentistry. everStick fiber reinforcements are made of silanated glass fibers in thermoplastic polymer and light curing resin matrix.

everStick products address the advantages of minimally invasive dentistry where the patient's own healthy tooth tissue is saved for as long as clinically possible. This also means that other treatment options remain available should the patient ever need them in the future.

## IPN - The heart of everStick® fibers

Proper bonding between the fibers and composite is the key factor for a successful treatment. Only everStick products have a unique, patented interpenetrating polymer network structure (IPN). Clinically this leads to superior bonding enabling reliable surface retained applications and perfect handling properties.

The significance of the IPN structure is that surfaces can be reactivated even after the final polymerisation. Reactivation is crucial for superior bonding when

- laboratory-manufactured restorations are cemented to teeth
- Fiber reinforced composite (FRC) devices are remodelled or repaired

The IPN structure makes the everStick products fundamentally different from any other fiber or composite materials available.

## Advantages of everStick® fibers

- Minimally invasive and reversible; leave other treatment options available
- Superior mechanical properties
- Unique patented bonding
- As strong as metal
- Elasticity close to that of dentin
- Extensive research data
- More than 300 publications and several dissertations
- Long term clinical data with excellent success ratio

# everStick® Product Range

## everStick®C&B

Content SKU  
 2 x 12 cm 900833  
 1 x 8 cm 900834  
 No. of fibers: 4000



## everStick®POST

Content SKU  
 Diameter 1.0, 10 posts 900828  
 Diameter 1.2, 10 posts 900829  
 Diameter 1.5, 10 posts 900830



## everStick®PERIO

Content SKU  
 2 x 12 cm 900822  
 1 x 8 cm 900832  
 No. of fibers: 2000



## everStick®NET

Content SKU  
 30 cm<sup>2</sup> 900818



## everStick®ORTHO

Content SKU  
 2 x 12 cm 900831  
 No. of fibers: 1600



## Stick®RESIN, 5 mL

SKU 900823



## G-aenial™ Universal Flo (2 mL/3.4 g per syringe)

Refill: 1 syringe, 20 Dispensing Tip III Plastic

Content SKU  
 Shade A2 004203  
 Shade A3 004204



## Stick®STEPPER

SKU 900825



## Stick®CARRIER

SKU 900826



## everStick®INTRO

SKU 900835  
 everStickC&B 8 cm  
 everStickPERIO 8 cm  
 StickREFIX D 1 silicone instr.



## everStick®COMBI

SKU 900836  
 everStickC&B 8 cm  
 everStickPERIO 8 cm  
 everStickNET 30 cm<sup>2</sup>  
 everStickPOST 1.2 5 posts  
 StickSTEPPER 1 hand instr.  
 Stick®REFIX D 1 silicone instr.



## everStick®POST INTRO

SKU 900827  
 everStickPOST 0.9 5 posts  
 everStickPOST 1.2 5 posts  
 everStickPOST 1.5 5 posts  
 StickRESIN 5 mL



## everStick®STARTER KIT

SKU 900819  
 everStickC&B 8cm  
 everStickPERIO 8cm  
 everStickNET 30cm<sup>2</sup>  
 everStickPOST Ø 1.2 5 posts  
 StickRESIN bottle 5 mL  
 G-aenial Universal Flo syringe 2 mL  
 1x StickREFIX D 1 silicone instr.  
 StickSTEPPER hand instr.  
 StickCARRIER hand instr.

