
**S1-E / S1-A**

**Sendoline®**

Operating Instructions  
Bruksanvisning

English / Svenska / Norsk / Dansk

Bruksanvisning  
Brugsanvisning

## Annex to IFU S1-E / S1-A

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EN - Annex to Operating Instructions S1-E / S1-A | Now you can also use our S1 Plus instruments with the S1-E / S1-A contra-angle handpiece!

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
Operating Instructions

# Content

<b>1</b>	<b>Before you begin</b> .....	<b>6</b>
1.1	Structure of the document .....	7
1.2	Service life .....	10
<b>2</b>	<b>Safety information</b> .....	<b>11</b>
<b>3</b>	<b>Technical description</b> .....	<b>13</b>
3.1	Task .....	13
3.2	Set up .....	13
3.3	Overview of files .....	14
3.4	Technical data .....	16
<b>4</b>	<b>Preparation</b> .....	<b>20</b>
4.1	Initial start-up and longer breaks in use .....	20

<b>5</b>	<b>Operation .....</b>	<b>21</b>
5.1	Replacing the instrument.....	21
5.2	Attaching and removing a file .....	22
5.3	Attaching the APEX clip .....	24
5.4	Connecting the instrument to an endometer .....	25
5.5	Recommendation for the treatment procedure.....	26
5.6	Mounting and removing the insulating sleeve .....	29
<b>6</b>	<b>Follow-up .....</b>	<b>31</b>
6.1	After each treatment session.....	31
6.2	At the end of the work day.....	32
<b>7</b>	<b>Conditioning .....</b>	<b>33</b>
7.1	Conduct pre-disinfection.....	33

7.2	Automated cleaning and disinfecting .....	34
7.3	Manual cleaning and disinfection .....	37
7.4	Manual maintenance .....	39
7.5	Sterilizing .....	43
<b>8</b>	<b>Spare parts and consumables .....</b>	<b>45</b>
<b>9</b>	<b>Storage and transport conditions .....</b>	<b>46</b>
<b>10</b>	<b>Disposal .....</b>	<b>48</b>

## 1 Before you begin ...

S1-E / S1-A complies state-of-the-art technical regulations. S1-E / S1-A complies with the ISO 14457 standard.

1. Read the operating instructions prior to using S1-E / S1-A.
2. Only use S1-E / S1-A for the applications described in the operating instructions.
3. Observe the applicable hygiene standards, occupational safety regulations and accident prevention measures for S1-E / S1-A.

S1-E / S1-A is used to hold and drive burr instruments for the purposes of rotary processing and is intended for the following uses in dentistry:

- Endodontics
- Endometry

None

This product is intended only for use by trained dental personnel in dental practices and laboratories.

**Intended Use**

**Contraindications**

**Target group**

## Warnings

## Instructions for use

### 1.1 Structure of the document

#### 1.1.1 Labeling of information

- To prevent injuries, please observe all warnings.

Warnings are labeled as follows:

 **DANGER!** indicates a danger **leading** to death or serious injury if not avoided.

 **WARNING!** indicates a danger that **may lead** to death or serious injury if not avoided.

 **CAUTION!** indicates a danger that **may lead** to injury if not avoided.

- To prevent material damage and additional expenses, please observe all instructions for use.

Instructions for use are labeled as follows:

**NOTICE!** indicates measures for the prevention of material damage.

**IMPORTANT:** indicates information on the avoidance of additional expenses and other important information.

**Tip:** indicates information for facilitating work.

### 1.1.2 **Formats and symbols**

The formats and symbols used in this document have the following meaning:

<ul style="list-style-type: none"><li>✓ Prerequisite</li><li>1. First action step</li><li>2. Second action step</li></ul> or <ul style="list-style-type: none"><li>➤ Alternative action</li></ul> <ul style="list-style-type: none"><li>↪ Result</li><li>➤ Individual action step</li></ul>	Requests you to do something.
Use of formats and symbols [ → 8].	Identifies a reference to another text passage and specifies its page number.
<ul style="list-style-type: none"><li>• List</li></ul>	Identifies a list.

## 1.2 Service life

When used as intended:

- Non-moving parts have a typical service life of 5 years
- Moving parts have a typical service life of 2 years

No warranty claim can be inferred here, as wear may occur earlier or later than indicated above depending on use, frequency of sterilization, and frequency of maintenance.

## Obligations of the user

### Preventing the spread of infections and cross contamination

### Instrument head/body overheating

## 2 Safety information

- Use only fault-free materials that do **not** deviate from the specified data [ → 16].
- Protect yourself, patients, and others against any foreseeable dangers. To do this, follow the safety information.
- Comply with the Intended use of the equipment.
- You should always keep these operating instructions within reach for further reference.

Prevent the spread of infections and cross contamination between patients, users, and third parties. Sterilize equipment after each patient. Take the appropriate hygiene measures, e.g. wear protective gloves. If the instrument is defective, the area around the instrument head or the instrument body may heat up, thus creating a risk of burning the hands or the patient's oral mucosa.

## **Malfunction or damage**

Discontinue use immediately in case of malfunction, unusual or different sounds or damage. Damaged instruments may cause injury. Notify the dental dealer or the manufacturer.

## **Repair**

Do **not** repair the instrument yourself.

## **Spare and accessory parts**

Use only original parts produced by the manufacturer.


If you have any questions, please contact your dental dealer or the manufacturer.

## 3 Technical description


### 3.1 Task

The instrument serves for transmitting the driving power of the electric motor or air motor (handpiece coupling according to ISO 3964) to the preparation tool and converts a rotational movement into a rotating movement with alternating direction.

### 3.2 Set up

A	Instrument head
B	Push button
	Opening of chuck system
D	Joint for clip

## Accessories

A	Insulating sleeve
B	Clip

### 3.3 Overview of files

Following files are to use for the application with S1-E / S1-A:

Designation	ISO/Taper
Sendoline – S1 File Standard 25/06	25/6
Sendoline – S1 File 15/06	15/6

<b>Designation</b>	<b>ISO/Taper</b>
Sendoline - S5 Rotary File #1	30/8
Sendoline - S5 Rotary File #2	30/6
Sendoline - S5 Rotary File #3	30/4
Sendoline - S5 Rotary File #4	25/4
Sendoline - S5 Rotary File #5	20/4
Sendoline - Big Apical File Green	35/4
Sendoline - Big Apical File Black	40/4
Sendoline - Big Apical File White	45/4
Sendoline - Big Apical File Yellow	50/4

Designation	ISO/Taper
Sendoline - S5 Retreatment File #1	30/12
Sendoline - S5 Retreatment File #2	30/8
Sendoline - S5 Retreatment File #3	30/6
Sendoline - S5 Retreatment File #4	30/4

### 3.4 Technical data

	S1-E	S1-A
Gear ratio	~ 150:1	~ 82:1
Maximum motor speed in rpm	~ 40000	~ 22000
Mean speed at maximum motor speed in rpm	~ 270	~ 270

	S1-E	S1-A
Motion	forward and backward rotating	forward and backward rotating
Clockwise rotation angle	~ 180°	~ 180°
Counterclockwise rotation angle	~ 30°	~ 30°
Periodicity of the motion in Hz	~ 10.7	~ 10.7
Clamping system	CA	CA
Internal cooling media	-	-
Light function	-	-
Apex measurement function	x	x

	S1-E	S1-A
Clips available for standard endometer connections	∅ 2mm	∅ 2mm
Instrument coupling	ISO 3964	ISO 3964

### File

	S1-E / S1-A
Shank diameter in mm	2.334 - 2.35
Maximum total length in mm	42
Maximum working diameter in mm (ISO 2157)	2.1

	<b>S1-E / S1-A</b>
Gripping length in mm	$\geq 11$
Standard	ISO 1797-1

This product bears the CE mark in accordance with the provisions of the Council Directive 93/42/EEC of June 14, 1993 concerning medical devices (MDD).

## **4 Preparation**

### **4.1 Initial start-up and longer breaks in use**

- Sterilize the instrument and accessories prior to startup.
- Clean and maintain the instrument after longer breaks in use.

## 5 Operation

**NOTICE!** Use only faultless files in order to prevent fatigue fractures of the files as far as possible.


**⚠ CAUTION!** Do not pull the patient's cheek back with the contra-angle handpiece while the motor is running. This would actuate the pushbutton, thus creating a risk of burning the patient's oral mucosa.

**⚠ CAUTION!** Use the S1-A only with an air-driven motor at a maximum speed of 22,000 rpm.

**IMPORTANT:** For application details and operating data of the different preparation tools, please refer to the information provided by the manufacturer.

### 5.1 Replacing the instrument

**⚠ CAUTION!** The instrument should only be fitted or removed when the motor is at standstill.

### Attaching the S1-E / S1-A

- ✓ The motor is at a standstill.
- Insert the instrument until it snaps into place.

### Removing the S1-E / S1-A

- ✓ The motor is at a standstill.
- Detach the instrument. Do not pull on the supply hose while doing this.


## 5.2 Attaching and removing a file

**NOTICE!** Use only faultless files in order to prevent fatigue fractures of the files as far as possible.

**NOTICE!** Only use files, which are approved by Sendoline.

**IMPORTANT:** The file is designed for reciprocal use.

### Attaching a file


- ✓ The drive has come to a stop.
- 1. Insert the file by gently rotating it until it snaps into place. In doing so, do **not** press the pushbutton.
- 2. Then pull and turn the file to check that it is firmly seated.


### Removing a file

- ✓ The file has come to a stop.
- Remove it while pressing and holding down the pushbutton.

### 5.3 Attaching the APEX clip

1. Attach the APEX clip to S1-E / S1-A.
2. Connect an endometer [ → 25].

## 5.4 Connecting the instrument to an endometer

The S1-E / S1-A makes it possible to connect a measuring cable for an endometer by means of a clip in order to check the exact position of the file in the root canal.

**NOTICE!** Use only endometers, which are approved in the respective country.

1. Please observe the operating instructions pertaining to the endometer being used.
2. Attach an APEX clip [ → 24].
3. Connect the measuring cable of the endometer to the clip's contact pin.
4. Perform the measurement in accordance with the equipment manufacturer's specifications.

Based on in vitro studies, we recommend defining the maximum preparation depth at the coronal end of the apex interval on the display of the endometer used.


### **Calibration**

- ✓ The manufacturer of the endometer stipulates that a calibration must be performed.
- 1. Insert a file into the contra-angle handpiece [ → 22].
- 2. Connect the cable of the file clamp (without the actual file clamp) onto the clip of the instrument.
- 3. Perform the calibration in accordance with the equipment manufacturer's specifications.

## **5.5 Recommendation for the treatment procedure**

**NOTICE!** Use only faultless files in order to prevent fatigue fractures of the files as far as possible.

- The different file systems require different working techniques. For this reason, please always observe the information provided by the file manufacturer.
- Work without exerting any pressure on the file.
- Prior to the first treatment, we recommend performing conditioning tests on resin blocks or extracted teeth.
- Due to the special characteristics of the nickel-titanium material, the crown down technique should be preferred when preparing the root canal.
- The ApexLocator function is used to support the root canal preparation process in combination with a treatment center with the apex measurement function. The number of X-rays can be reduced with S1-E / S1-A. However, you must always take at least one X-ray to determine the depth for preparatory work. A clinical evaluation including knowledge of the anatomy of root canals is important for interpreting the results.

- When performing the APEX measurement, wear appropriate insulated gloves. We recommend performing the treatment with a cofferdam, while using an insulating sleeve. This will prevent inaccurate measurements caused by undesired leakage currents. During the measurement, the instrument must not come into contact with the patient's mucosa, metallic tooth restorations, or the mucosal electrode.
- If the limit torque is exceeded, a gentle creaking may be heard and vibrations felt. Remove the file by gently pulling it toward the coronal end. If this is not possible, switch the drive to counterclockwise rotation and pull the file out of the channel.  
**IMPORTANT:** Check the file afterwards for damage or deformation (untwisting) and replace it, where necessary.

## 5.6 Mounting and removing the insulating sleeve


**NOTICE!** Insulating sleeves are only used once. Change the insulating

sleeve after each patient.

### Fitting an insulating sleeve

- ✓ You have sterilized the insulating sleeve.
  - ✓ The motor is at a standstill.
  - ✓ A file is **not** clamped in place.
1. Twist the insulating sleeve carefully this way and that and pull the insulating sleeve fully over the instrument.
  2. Attach the file [ → 22].

### **Removing an insulating sleeve**

✓ The motor is at a standstill.

1. Remove the file.
2. Remove the insulating sleeve.
3. Dispose of the insulating sleeve [ → 48].

## 6 Follow-up

### 6.1 After each treatment session

**NOTICE!** Condition immediately, or at the latest, one hour after treatment.

**NOTICE!** Only use W&H Service Oil F1, MD-400.

✓ Wear appropriate protective clothing.

1. Remove the file with tweezers.
2. Pre-disinfect directly at the treatment center [ → 33].
3. Detach the instrument from the drive.
4. Transport the instrument to the hygiene room in a suitable transport container.
5. Conduct automatic conditioning [ → 34]. Manual conditioning [ → 37] is possible in exceptional cases if the national/local regulations are followed.

6. Apply spray to the instrument [ → 39].
7. Sterilize the instrument and accessories [ → 43].

## 6.2 At the end of the work day

- Apply spray to the instrument [ → 39].

**NOTICE!** Do not leave any instruments on the motor overnight, in order to prevent oil from leaking into the electric motor. Never lubricate the electric motor.

## 7 Conditioning


### 7.1 Conduct pre-disinfection

- ✓ Wear appropriate protective clothing.
- ✓ All disinfectants must be approved in your country and have proven bactericidal, fungicidal and virucidal properties. Only use disinfectants with **no** protein-fixing properties.
  1. Spray the surface with disinfectant.
  2. Wipe the disinfectant off with a cloth.
- ✎ For further conditioning, the instrument should be dry and free of residue.

Please observe the manufacturer's instructions for using instrument disinfectants.

## 7.2 Automated cleaning and disinfecting ...

### 7.2.1 ... with a DAC UNIVERSAL

We recommend using **Sirona DAC UNIVERSAL** for automated cleaning, disinfection, care, and unwrapped sterilization.

For further details, refer to the operating instructions supplied with the unit.

- ✓ The pushbutton chuck is manually maintained.
  - ✓ The instrument is conditioned using the DAC UNIVERSAL.
1. Check whether the instrument is clean after conditioning under good lighting (min. 500 lux) and color rendering index (min. 80 Ra).
  2. If it is still dirty, repeat the process.
    - ↳ For further conditioning, the instrument should be dry and free of residue.

3. Pack the instrument in packaging material suitable for sterilization and storage. e.g. paper/plastic composite packaging.
4. Sterilize the instrument [ → 43].

### 7.2.2 ... with cleaning and disinfection equipment

The instrument can also be cleaned and disinfected in a suitable piece of cleaning and disinfection equipment.

The cleaning and disinfection equipment used must be approved by its manufacturer for the cleaning and disinfection of dental instruments and comply with EN ISO 15883-1 (e.g. 95°C (203°F) and 10 min. holding time).

For further details, refer to the operating instructions supplied with the unit.

- ✓ The instrument is conditioned with a cleaning and disinfection device.

1. Check whether the instrument is clean after conditioning under good lighting (min. 500 lux) and color rendering index (min. 80 Ra).
2. If it is still dirty, repeat the process.
  - ↳ For further conditioning, the instrument should be dry and free of residue.
3. Blow the instrument out with max. 3 bar.
4. Manual care of mechanical parts [ → 39].
5. Manual care of the push button chuck.
6. Pack the instrument in packaging material suitable for sterilization and storage. e.g. paper/plastic composite packaging.
7. Sterilize the instrument [ → 43].

### 7.3 Manual cleaning and disinfection

**IMPORTANT:** Manual conditioning is possible in exceptional cases if the national/regional regulations are followed. The national/regional regulations are to be checked before.


**NOTICE!** Condition immediately, or at the latest, one hour after treatment.

**NOTICE! Never** clean in an ultrasound bath!

**NOTICE! Never** immerse in disinfectant solution!

**IMPORTANT:** Use a soft, clean, and disinfected brush for cleaning.

- ✓ Wear appropriate protective clothing.
- ✓ All disinfectants must be approved in your country and have proven bactericidal, fungicidal and virucidal properties. Only use disinfectants with **no** protein-fixing properties.

1. Brush the instrument under running water (< 38 °C, < 100 °F, at least drinking water quality) and good lighting (min. 500 lux) and color rendering index (min. 80 Ra) until no more dirt can be seen, for at least 10 seconds.
2. Flush the drive channels with spray.
3. Conduct thermal disinfection or unwrapped steam sterilization.
4. Maintain mechanical parts manually [ → 39].
5. Maintain the push button chuck manually.
6. Pack the instrument in packaging material suitable for sterilization and storage. e.g. paper/plastic composite packaging.
7. Sterilize the adapter [ → 43].

## 7.4 Manual maintenance

### 7.4.1 Maintenance of mechanical parts

#### Intervals


- At least every noon and evening
- Prior to each sterilization
- After every thermal disinfection without integrated maintenance

**NOTICE!** Only use W&H Service Oil F1, MD-400.

#### Process

##### With W&H Service Oil F1, MD-400

- ✓ The spray adapter is disinfected.
1. Fit the spray cap with the nozzle onto the spray can.
  2. Shake the spray can before each application.

3. Insert the nozzle into the Instrument.
4. Hold the instrument firmly.
5. Spray the instrument for about 1 second.  
**IMPORTANT:** Hold the spray can upright.
6. Wipe any spray that comes out with a disinfection cloth.
7. Repeat the process until the spray is clear.


#### **With W&H Assistina**

- Observe the operating instructions of the manufacturer.

### 7.4.2 Care of the push button chuck


#### **Interval**

- At least once a week

**Only with W&H Service Oil F1, MD-400**

✓ The instrument is clean and disinfected.


1. Fit the spray cap on the spray can.
2. Hold the instrument firmly.
3. Press the tip of the spray nozzle firmly into the chuck system.
4. Spray into the instrument for about 1 second.  
**IMPORTANT:** Hold the spray can upright.
5. Wipe any spray that comes out with a disinfection cloth.

### Only with W&H Assistina

✓ The instrument is clean and disinfected.

1. Fit the adapter onto the maintenance attachment.
2. Direct the spray exit nozzle of the adapter downwards.
3. Press the bur side of the instrument onto the adapter from below.
4. Start the Assistina.
5. Press the instrument onto the adaptor for about 10 seconds.
6. Remove the instrument from the adapter.
7. Wipe any spray that comes out with a disinfection cloth.
8. Close the cover of the Assistina and allow it to run for the remainder of the maintenance cycle (approx. 25 seconds).

## 7.5 Sterilizing

- ✓ The instrument is clean and disinfected.
- ✓ Instrument maintenance is complete.
- ✓ The instrument can be sterilized in packaging suitable for sterilization and storage: paper/plastic composite packaging according to ISO 11607.
- Sterilize the instrument in the steam sterilizer with saturated water vapor.

Overpressure: 2.04 bar (29.59 psi)

Temperature: 134° C (274° F)

Holding time: 3 min

Steam sterilizers that meet the requirements of EN 13060, class B or S, and are also suitable for the sterilization of straight/contra-angle handpieces are approved.

**NOTICE!** Do not exceed 140°C (284°F), even during the drying phase.

#### **After sterilizing**

1. Remove the instrument from the steam sterilizer immediately.  
**⚠ CAUTION!** The instrument is hot. Risk of burns!  
**NOTICE!** Do **not** attempt to accelerate the cooling process by immersing the instrument in cold water. This can damage your instrument.
2. Store all instruments so that they are protected from contamination.
3. Sterilize again once the storage period has elapsed.

## 8 Spare parts and consumables

	REF			REF
Service Oil F1, MD-400 (package of 6 oil spray cans, per 400 ml)	10940021		S1 apex clip Ø 2.0 mm	65 83 566
Spray cap with nozzle	02038200		S1 insulating sleeve (for S1-E/S1-A)	65 83 541
W&H spray cap with spray nozzle for chuck system	02036100		W&H Assistina adapter for chuck system	02693000

## 9 Storage and transport conditions


Protect from  
moisture


Sensitive  
contents


Relative  
humidity

Air pressure

Temperature

After a severe change in  
temperature, allow sufficient  
time for acclimation.

## 10 Disposal

- According to current information, the product does not contain any substances that are hazardous to the environment.
- Disinfect the product prior to disposal.
- Observe the applicable disposal regulations for your area.

We reserve the right to make any alterations which may be required due to technical improvements.

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Beteckning	ISO/Taper
Sendoline - S1 File Standard 25/06	25/6
Sendoline - S1 File Small 15/06	15/6
Sendoline - S5 Rotary File #1	30/8
Sendoline - S5 Rotary File #2	30/6
Sendoline - S5 Rotary File #3	30/4
Sendoline - S5 Rotary File #4	25/4
Sendoline - S5 Rotary File #5	20/4
Sendoline - Big Apical File Green	35/4

<b>Beteckning</b>	<b>ISO/Taper</b>
Sendoline - Big Apical File Black	40/4
Sendoline - Big Apical File White	45/4
Sendoline - Big Apical File Yellow	50/4
Sendoline - S5 Retreatment File #1	30/12
Sendoline - S5 Retreatment File #2	30/8
Sendoline - S5 Retreatment File #3	30/6
Sendoline - S5 Retreatment File #4	30/4

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Betegnelse	ISO/taper
Sendoline - S1 File Standard 25/06	25/6
Sendoline - S1 File Small 15/06	15/6
Sendoline - S5 Rotary File #1	30/8
Sendoline - S5 Rotary File #2	30/6
Sendoline - S5 Rotary File #3	30/4
Sendoline - S5 Rotary File #4	25/4
Sendoline - S5 Rotary File #5	20/4
Sendoline - Big Apical File Green	35/4

<b>Betegnelse</b>	<b>ISO/taper</b>
Sendoline - Big Apical File Black	40/4
Sendoline - Big Apical File White	45/4
Sendoline - Big Apical File Yellow	50/4
Sendoline - S5 Retreatment File #1	30/12
Sendoline - S5 Retreatment File #2	30/8
Sendoline - S5 Retreatment File #3	30/6
Sendoline - S5 Retreatment File #4	30/4

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Betegnelse	ISO/Taper
Sendoline - S1 File Standard 25/06	25/6
Sendoline - S1 File Small 15/06	15/6
Sendoline - S5 Rotary File #1	30/8
Sendoline - S5 Rotary File #2	30/6
Sendoline - S5 Rotary File #3	30/4
Sendoline - S5 Rotary File #4	25/4
Sendoline - S5 Rotary File #5	20/4
Sendoline - Big Apical File Green	35/4

<b>Betegnelse</b>	<b>ISO/Taper</b>
Sendoline - Big Apical File Black	40/4
Sendoline - Big Apical File White	45/4
Sendoline - Big Apical File Yellow	50/4
Sendoline - S5 Retreatment File #1	30/12
Sendoline - S5 Retreatment File #2	30/8
Sendoline - S5 Retreatment File #3	30/6
Sendoline - S5 Retreatment File #4	30/4

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**Sendoline AB**