CHARISMA® flow

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Made in Germany

Distributed in USA/Canada exclusively by:
Heraeus Kulzer, LLC
300 Henkle Way,
South Bend, IN 46614-2517
1-800-431-1785

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CHARISMA® flow

Charisma flow is a flowable, light curing, radiopaque hybrid composite used for adhesive, tooth-coloured anterior and posterior restorations. Due to its low viscosity Charisma flow provides excellent coverage of the tooth structure.

Composition: Charisma Flow is based on multifunctional methacrylate monomers (EBADMA/TEGDMA) and contains approximately 62% by weight or 38% by volume inorganic fillers such as Ba-Al-F silicate glass and SiO₂. These particles size is between 0.005µm and 5µm.

Shades: A1, A2, A3, A3.5, B2, B3, OA2, OA4, OB2, SB, Baseline

Applications:

- Fissure sealing
- Enlarged fissure sealing
- Class V fillings
- Minimally invasive Class I and II fillings in areas not subjected to masticatory forces
- Minimally invasive Class II fillings
- Small contour and shade adjustment on the enamel and dentine
- Minimally invasive Class II fillings
- Restorations peu étendue de restaurations directes et indirectes en combinaison avec un adhésif adéquat
- Contentions dentaires
- Revêtement de cavités
- Collage de brackets

Kind of application:

- Clean the tooth structure for treatment by cleaning it with a fluoride free polishing paste. Before treatment, while the teeth are still moist, select the Charisma flow shade using the Charisma shade guide.
- Remove any caries and complete minimally invasive preparation of the tooth structure according to the general guidelines for adhesive fillings. Then clean and dry the cavity. The use of a rubber dam is recommended.
- Apply an adhesive system (e.g. Bond® Self Etch, GLUMA® Comfort® Bond + Desensizer, GLUMA® Solid Bond) should be applied to the enamel and dentin in accordance with the manufacturer's instructions before inserting the filling material. Avoid any contamination of conditioned or adhesive-treated surfaces with moisture, blood or saliva.
- Apply Charisma flow in thin layers (max. 2mm, Baseline max. 1 mm) to the cavity. When applying Charisma flow ensure that it flows slowly and evenly to prevent air bubbles – apply a gentle stream of air to facilitate maximal marginal contact. If necessary, remove any air bubbles with a probe before polymerisation.
- Apply each layer according to the table below using a suitable light-curing unit (e.g. Translux® Power Blue®) for 20 or 40sec. A dispersion layer forms during polymerisation. This layer should not be touched or removed if other composite layers are to be applied.

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