## CONTENT

<b>EN (English)</b> Instruction Manual GenTek™ TiBase.....	2
<b>BG (Български)</b> Ръководство за употреба База GenTek™ TiBase.....	5
<b>CS (Česky)</b> Návod k použití GenTek™ TiBase.....	8
<b>DA (DANISH)</b> Brugervejledning GenTek™ TiBase.....	11
<b>DE (Deutsch)</b> Gebrauchsanweisung GenTek™ TiBase.....	14
<b>EL (Ελληνικά)</b> Εγχειρίδιο οδηγιών Στήριγμα τιτανίου GenTek™.....	17
<b>ES (Español)</b> Manual de instrucciones GenTek™ TiBase.....	20
<b>ET (Eesti)</b> Kasutusjuhend GenTek™ TiBase.....	23
<b>FI (Suomi)</b> Käyttöohje GenTek™ TiBase -kanta.....	26
<b>FR (Français)</b> Mode d'emploi Base titane GenTek™.....	29
<b>HR (Hrvatski)</b> Priručnik GenTek™ titanska baza.....	32
<b>HU (Magyar)</b> Kezelési útmutató GenTek™ TiBázis.....	35
<b>IT (Italiano)</b> Istruzioni d'uso GenTek™ TiBase.....	38
<b>LT (Lietuvių k)</b> Naudojimo instrukcija GenTek™ TiBase.....	41
<b>LV (Latviski)</b> Lietošanas rokasgrāmata GenTek™ TiBase.....	44
<b>NL (Nederlands)</b> Gebruikershandleiding GenTek™ TiBase.....	47
<b>NO (Norsk)</b> Instruksjonshåndbok GenTek™ TiBase.....	50
<b>PL (Polski)</b> Instrukcja stosowania TiBase GenTek™.....	53
<b>PT (Português)</b> Instruções de uso Base de titânio GenTek™.....	56
<b>RO (Română)</b> Manual de instrucțiuni Bază din titan GenTek™.....	59
<b>SK (Slovak/Slovakian)</b> Návod na použitie Titánová báza GenTek™.....	62
<b>SL (Slovenščina)</b> Navodila za uporabo GenTek™ TiBase.....	65
<b>SV (Svenska)</b> Bruksanvisning GenTek™ TiBase.....	68
<b>TR (Türkçe)</b> Kullanım Kılavuzu GenTek™ TiBase.....	71

## Indications

### Intended Use / Intended Purpose

The GenTek™ TiBase along with the final prosthetic screw is intended to support screw and cement-retained single-unit, multi-unit, and full-arch prosthetic restorations. It is used as an interface between a dental implant and an abutment restoration.

The gluing screw is intended for laboratory use during the cementing process. The gluing screw is used to block the screw channel from the potential inflow of glue while the prosthetic restoration is being glued to the TiBase, prior to final placement of the TiBase and screw into the mouth.

### Device description

Each GenTek™ TiBase contains a straight screw channel and is packaged with the final prosthetic screw and a black gluing screw. The final prosthetic screw for each system is also available as an individually packaged component. GenTek™ TiBase, prosthetic screw, and gluing screw are for single use only. GenTek™ TiBase and prosthetic screw are provided non-sterile and intended to be cleaned and sterilized prior to placement into the patient's mouth.

The GenTek™ portfolio is comprised of implant-level and abutment-level tibases.

#### Implant-Level

The Implant-Level GenTek™ Tibases contain an overall post height of 4.7mm with a horizontal cutting groove to indicate a height adjustment of 4.0mm. Additionally, the GenTek™ Tibases are Sirona-compatible and contain an anti-rotation notch designed to prevent the restoration from rotating.

All Implant-Level GenTek™ Tibases are available with genuine connections for the following Zimmer Biomet Dental implant systems: Certain®, External Hex, TSV® / Trabecular Metal™, and Ezetic®.

#### Abutment-Level

Abutment-Level GenTek™ Tibases are available for the Zimmer Biomet Dental Low Profile and Tapered Abutment Systems.

For specific dental implant compatibility, please refer to the information on the individual product label and the GenTek™ components catalog. All GenTek™ Tibases are compatible with Zimmer Biomet Dental prosthetic tooling as outlined in the GenTek™ Components Catalog – please refer to the components catalog located on the Zfx website.

### Materials

GenTek™ Tibases are manufactured from biocompatible titanium alloy grade 5 ELI (TiAl6V4 ELI).

## Procedure

### × Designing the Custom Restoration

The GenTek™ TiBase can be processed via a digital or traditional workflow utilizing GenTek™ Scan Bodies. The custom restoration must be constructed using dental design software. The geometry of the TiBase abutment is contained in the library file. Minimum wall thickness requirements should be controlled according to the material supplier's recommendations.

### × Milling a Custom Restoration

The custom restoration must be manufactured with milling equipment adequate for the designed geometry and material of the restoration.

### × Post-Processing the Custom Restoration

Follow the post-processing steps in accordance with the instructions issued by the material supplier.

### × Preparing the Custom Restoration for Cementing

Once the custom restoration is ready for cementing, refer to the material supplier's instructions for use for proper cementing process.

### × Preparing the GenTek™ TiBase for Cementing

Before cementing the restoration (crown or framework) to the abutment, the gluing surface of the TiBase should be thoroughly cleaned as outlined in the **Cleaning and Sterilization** section of this IFU. Blast GenTek™ TiBase in accordance with the instructions issued by the cement supplier.

#### GenTek™ Tibases (Implant-Level)

The post of the GenTek™ Implant-Level Tibases contain a groove which indicates a 4 mm post height. These abutments should not be cut back below this groove. The maximum post height of the GenTek™ Implant-Level Tibases is 4.75 mm and contains a cutting groove which indicates a 4 mm post height. The abutment can be left at its maximum post height of 4.75 mm or it can be cut at the 4.0 mm groove that best fits the case so that the crown can be digitally designed from the corresponding library file.

#### GenTek™ Tibases (Abutment-Level)

The GenTek™ Abutment-Level Tibases do not contain a reference height via a groove and therefore should not be cut.

### × Cementing

Follow cementing procedures in accordance with the instructions issued by the cement supplier. In order to cement the restoration to the GenTek™ TiBase, insert the base into an implant analog and fix it in place with the gluing screw (finger-tighten). The gluing screw should be kept in place throughout the cementing process. Apply the cement in a circular fashion to the top edge of the contact surface of the GenTek™ TiBase and to the internal geometry of the restoration. Rotate the restoration gently as you push it over the GenTek™ TiBase until resistance is felt. For the Implant- and Abutment-Level tibases, the rotation during assembly of the restoration does not apply. Turn the restoration until it snaps into the notch feature on the TiBase, and then proceed to bring the restoration into its final position.

**Note:** Avoid raising the restoration after the final position has been reached, as this can compromise the integrity of the bond. The custom restoration must fit the GenTek™ TiBase without gaps and voids. Any traces of cement should be removed immediately. The restoration must be held in position during cement removal.

## × Polishing

Once the cement has set, remove the excess cement with rubber polishing tools. When polishing, fix the TiBase/restoration in an implant analog to prevent damaging the connection geometry.

## × Placing / Removing the TiBase Restoration

### Placing in a Laboratory Setting

When placing the TiBase restoration onto the analog in the model, use the corresponding screw to affix the restoration to the analog and finger-tighten using the corresponding driver – refer to Table 1 for driver compatibility by connection.

### Placing in a Clinical Setting

1. Remove any provisional abutment, cover screw, healing abutment or similar from the implant upon which the TiBase is intended to be placed.
2. Confirm that the selected TiBase is the same size as the implant or abutment-level connection.
3. Ensure that the seating surface of the implant or abutment-level connection is clean, free of residue, bone, or soft tissue and is dry when mounting.
4. Affix the Tibase restoration into the implant or abutment-level connection using the appropriate GenTek final prosthetic screw and Zimmer Biomet Dental driver. Refer to the below table for placement torque recommendations for the final prosthetic screw and for Zimmer Biomet driver compatibility:

**Note: Unlike the Certain® and External Hex TiBase systems, the TSV® and Eztetic® TiBase engaging systems form a friction-fit with the mating implant upon final assembly, requiring the use of additional instrumentation to remove the abutment from the implant.**

Connection	Final Prosthetic Screw Torque	Driver for Final Prosthetic Screw
<b>GenTek™ TiBase (Implant-Level)</b>		
Certain® (Gold-Tite® Screw, hex drive)	20 Ncm	PHD02N, PHD03N Tip: RASH3N, RASH8N
External Hex (Gold-Tite® Screw, square drive)	32 Ncm	PSQD0N, PSQD1N Tip: RASQ3N, RASQ8N
TSV®/Trabecular Metal™	30 Ncm	HXGR1.25, HXLGR1.25, HX1.25D, HXL1.25
Eztetic®	30 Ncm	HXGR1.25, HXLGR1.25, HX1.25D, HXL1.25
<b>GenTek™ TiBase (Abutment-Level)</b>		
Tapered	20 Ncm	HXGR1.25, HXLGR1.25, HX1.25D, HXL1.25
Low Profile	10 Ncm	PHD02N, PHD03N Tip: RASH3N, RASH8N

**Table 1**

\*For more information refer to the Zfx GenTek™ Product Catalog.

## × Removal of Friction-Fit Connections

Should the implant-level TiBase (TSV® / Trabecular Metal™ and Eztetic®) restoration require disassembly from the implant upon loosening and removal of the screw utilizing a GemLock Tool, an abutment removal tool may be required. Refer to GenTek™ Friction-Fit Abutment Removal Tool IFU-4097 for abutment removal instructions.

Removal tool for narrow platform (NP – Eztetic)

Removal tool for standard platform (SP – TSV/TM)

## Cleaning and Sterilization

The GenTek™ TiBase is non-sterile. Prior to placement of the GenTek™ TiBase / restoration and prosthetic screw into the implant, the GenTek™ TiBase assembly must first be cleaned and sterilized.

### Step 1: Manual Cleaning

1. Rinse the instruments under running cold tap water for 2 minutes.
2. Hand wash instruments in a bath with demineralized water for 5 minutes. Remove debris with a soft brush. Cavities should be rinsed with a syringe. This step should be repeated until all debris is removed.
3. Clean the instrument in an ultrasonic bath using a mildly alkaline cleaner (1.5%) at a frequency of at least 40 kHz for 20 minutes.
4. Rinse thoroughly with demineralized water for 1 min.

### Step 2: Sterilization

One of the following sterilization methods is recommended to be used:

#### × Method 1:

Full cycle pre-vacuum steam sterilization at a temperature of 134°C (273.2°F) for 3 minutes, 4 pulse minimum.

#### × Method 2: (alternative to method 1)

Full cycle gravity steam sterilization at 134°C (273.2°F) for 3 minutes minimum.

### Step 3: Drying

After sterilization treatment it is recommended to dry the products at a temperature of maximum 134°C (273.2°F) for 10 minutes. The products are ready for use 5 minutes after completion of the drying process.

**Note: Follow the loading instructions issued by the sterilizer manufacturer.**

## Storage and Handling

The GenTek™ TiBase should be stored at room temperature.

## Contraindications

The GenTek™ TiBase is not intended for restorations requiring an angle more than 20° to the axis of the implant and post height less than 4 mm.

## Known Side Effects

Allergies to the alloy or contents of the alloy may very rarely occur.

## Warnings

- × The GenTek™ TiBase must only be used and handled by dental professionals who have approved training in dentistry.
- × GenTek™ TiBase may only be used for Zimmer Biomet Dental implant systems.
- × The black gluing screw is not intended to be used in the patient's mouth.
- × The GenTek™ TiBase and prosthetic screw are small and, therefore, they must be handled with caution to avoid swallowing or inhaling by the patient.
- × The use of the prosthetic screw torque higher than the recommended value may result in damage to the GenTek™ TiBase, the prosthetic screw and/or the implant. The use of torque values lower than those recommended may result in loosening of the GenTek™ TiBase and/or the prosthetic screw.
- × The screwdriver must be in good condition. A worn screwdriver can lead to over-rotation of the hex in the screw head which can result in the screw not being fully tightened and/or not being able to be removed.
- × For any lab work and any try-in of the GenTek™ TiBase, a lab screw must be used. The included prosthetic screw may only be used to permanently fix the finished restoration.
- × Abutment level TiBases are not intended for single-unit restorations.
- × The gluing screw is not intended to be used when the screw channel is angulated.

## Explanation of symbols

Product is packaged non sterile

Do not use if package is damaged

Manufacturer Information: Address of legal manufacturer's facility

Don't use twice

Recommended torque to be applied for final placement of the included screw

Recommended to fix the included screw hand-tight

LOT-Number, Symbol followed by Production Date of the lot

Medical Device

Reference Number, Symbol followed by Manufacturer's Item Number

Quantity / Unit of package content

Unique Device Identifier

Refer to Instruction for use:  
Document is available for download on Zfx-Homepage

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Zfx GmbH