Instructions for use

Product description:

Ufi Gel SC is an A-silicone based permanently soft relining material. Its autopolymerising, cold-curing properties assure that **Ufi Gel SC** is simple and quick to use for direct as well as indirect relining. The set contains an adhesive liquid to achieve a firm bond between denture and relining material, the actual relining material, as well as glazing base and catalyst for smoothing surfaces after grinding. **Ufi Gel SC** is suitable for all PMMA-based dentures. A newly made **Ufi Gel SC** relining can be repaired and corrected.

Indications:

- Permanently soft relining for total and partial dentures
- to eliminate pressure spots
- to dam the palatal vibrating line
- in cases of flabby ridges and/or insufficient adhesion
- to cushion sharp-edged alveolar processes
- to support the healing process in implantology

Method of application:

1. Preparation and coating of the denture:
- Clean denture thoroughly with a brush and dry. Always remove any existing previous relinings. Grind down all areas to be relined by at least 1 mm to allow a later relining of at least 2 mm thickness. Important: chamfer preparation (shoulder preparation) of the rims of all areas to be relined. Prepare vestibular, lingual and dorsal rims in particular in the form of a chamfer. Round off the rims of the denture. Remove grinding dust with a dry brush or dry, oil-free air. Do not use water or a steam jet. Then degrease the denture surface by wetting a dust-free cloth with pure alcohol (min. 90 %) and wiping the surface clean. Then let aerate for 1 min. Humidity, additives (i. e. as in cavity cleaning agents), saliva, monomer-containing synthetic liquids or touching will diminish the adhesion of **Ufi Gel SC** to the denture.

Attention: Close bottle cap of the adhesive for the brush cap prior to first use.

Do not use the brush cap for the alcohol!

Apply the adhesive to the prepared surfaces which are to retain **Ufi Gel SC**. The surface has to be covered completely and evenly with the adhesive.

Attention: Close bottle immediately after use and store in an upright position.

Let aerate for 1 min.

Insert the **Ufi Gel SC** cartridge into a **VOCO** dispenser (type 2) or any dispensing gun with a suitable plunger. Slight differences in the filling levels of the chambers may occur for technical reasons. If this happens, remove the transport cap from the cartridge and extract material until both chambers deliver uniformly. When changing tips/caps, always ensure that the discharge openings of the **Ufi Gel SC** cartridge are clear. Then fix and lock the mixing tip with a 90° turn. The material is extruded and automatically mixed in the correct ratio by pressing the dispensing gun lever. After use, leave the mixing tip on the cartridge or use the transport cap.

Apply the mixed material first to the lingual/dorsal rim, then vestibular. Then coat the base. Make sure to cover all areas carrying adhesive with silicone. On areas to be relined, **Ufi Gel SC** should be at least 2 mm thick. **Ufi Gel SC** has a working time of 1 min 45 s from commencement of mixing.

2a. Direct relining in the mouth:

1 min after coating with **Ufi Gel SC**, the denture can be inserted into the patient's mouth.

Have the patient gently close the mouth in normal occlusion (1 min).

Then the patient should carry out masticatory and swallowing movements (5 min).

After 6 min in total, the denture can be removed from the mouth.

Important: Strictly observe the above-stated time intervals.

Finishing should be carried out 10 min after removal of the relined denture from the mouth.

Any earlier finishing should be carried out with extreme care.

2b. Indirect relining in the laboratory:

When relining is carried out in the laboratory, coat the denture and then isolate the dampened plaster model with a formaldehyde-free alginate-based separating agent for synthetic materials (not recommended: silicone-based separating agent) according to its instructions for use. Place the denture coated with **Ufi Gel SC** on the model and fix both in the relining device or casting flask. Shape functional rims with a metal spatula or a clean plaster knife. After shaping, let **Ufi Gel SC** set for 2 - 5 min, then place the relining device (or casting flask) loaded with the denture and the model in the pressure unit for 15 min at 40 °C - 45 °C. Afterwards leave the relining device to cool for 20 min.

Attention: The denture, plaster model or relining device may not be worked on with a steam jet in any phase.

3. Finishing the relining:

Any excess material can be removed with fine sharp (cuticle) scissors or a scalpel. Finishing of the transition area: If necessary, carry out grinding after 30 min, using polishing disks (REF 2049) or abrasive strips. Any earlier grinding process needs to be carried out with great caution. The **Ufi Gel SC** set contains 1 polishing disk REF 2049. Abrasive strips are available in varying grain sizes. The most suitable sizes are 80 (coarse) and 120 (medium). Use medium to low rotational speed (5,000 - 8,000 RPM) and apply only very light pressure when finishing the transition area. This prevents heat generation in the relining as well as excessive abrasion of the denture material.

Then clean with alcohol, allow any warmed silicone material to cool down and apply glazing.

4. Corrections:

Immediate correction of the relining is possible. Rinse the relining with water and dry, then apply fresh material, shape it and let it set in the mouth for 7 min or in the pressure unit for 5 min.

Any later correction (within two weeks) of the relining can be carried out after removal of the top layer of **Ufi Gel SC**.

5. Sealing:

Use glazing to seal and smoothen the surfaces treated in step 3. For better dosing, hold the bottle vertically with the dropper down and let air bubbles rise to the top. Mix an equal number of drops of glazing base and catalyst with a brush on a mixing slab until homogeneous. Apply the mixture thinly on the corresponding margins with a brush, i. e. on the fresh, dry cut and ground areas. The working time of **Ufi Gel SC** glazing is 1 min: 30 s. The setting time at room temperature is 10 min. Cooler conditions prolong the setting time. The setting time can be shortened to 3 min by using a pressure unit (40 °C - 45 °C).

Warnings, precautionary measures:

- To prepare the denture and improve the adhesion, additional roughening of the denture material carried out using sandblasting. Afterwards, wipe with pure alcohol.
- Process **Ufi Gel SC** at room temperature.
- Storage of **Ufi Gel SC** in a refrigerator prolongs the working and setting times.
- Always observe the setting times until final setting to achieve an optimal adhesive bond.
- Seal only the rough, finished areas (functional rim) with glazing to prevent accumulation of food residue or discoloration.
- Prevent any contamination with methacrylates (e.g. pressure units contaminated after denture fabrication and methacrylate-containing solvents for the removal of relining materials).
- Contamination may lead to incomplete curing of **Ufi Gel SC**/glazing and an irritating taste in the mouth.
- Close glazing bottles immediately after use:
 - Glazing base = white cap
 - Glazing catalyst = red cap
- For finishing, use the enclosed polishing disc, coarse, brown, REF 2049 or other suitable burs or polishing discs designed for flexible materials. Do not inhale silicone or synthetic dust when finishing the denture or the relining. Wearing a protective mask is recommended.
- **Ufi Gel SC** contains modified polydimethylsiloxanes and platinum catalyst. The adhesive contains butanone and methacrylates. **Ufi Gel SC** should not be used in case of known hypersensitivities (allergies) to any of these ingredients.
- Cleaning of the relined denture should be carried out with a special soft toothbrush or toothbrush and an abrasion-free cleanser (e. g. neutral soap). Rinse thoroughly afterwards.

Storage:

Store at temperatures between 4 °C - 23 °C.
Do not use after expiry date.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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3. Finishing der Relining:

Any excess material can be removed with fine sharp (cuticle) scissors or a scalpel.

Finishing of the transition area: If necessary, carry out grinding after 30 min, using polishing disks (REF 2049) or abrasive strips.

Any earlier grinding process needs to be carried out with great caution.

The **Ufi Gel SC** set contains 1 polishing disk REF 2049. Abrasive strips are available in varying grain sizes. The most suitable sizes are 80 (coarse) and 120 (medium). Use medium to low rotational speed (5,000 - 8,000 RPM) and apply only very light pressure when finishing the transition area. This prevents heat generation in the relining as well as excessive abrasion of the denture material.

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Use glazing to seal and smoothen the surfaces treated in step 3.

For better dosing, hold the bottle vertically with the dropper down and let air bubbles rise to the top. Mix an equal number of drops of glazing base and catalyst with a brush on a mixing slab until homogeneous.

Apply the mixture thinly on the corresponding margins with a brush, i. e. on the fresh, dry cut and ground areas. The working time of **Ufi Gel SC** glazing is 1 min: 30 s. The setting time at room temperature is 10 min. Cooler conditions prolong the setting time. The setting time can be shortened to 3 min by using a pressure unit (40 °C - 45 °C).

Any earlier finishing should be carried out with great caution.

The **Ufi Gel SC** set contains 1 polishing disk REF 2049. Abrasive strips are available in varying grain sizes. The most suitable sizes are 80 (coarse) and 120 (medium). Use medium to low rotational speed (5,000 - 8,000 RPM) and apply only very light pressure when finishing the transition area. This prevents heat generation in the relining as well as excessive abrasion of the denture material.

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5. Sealung:

Use glazing to seal and smoothen the surfaces treated in step 3.

For better dosing, hold the bottle vertically with the dropper down and let air bubbles rise to the top. Mix an equal number of drops of glazing base and catalyst with a brush on a mixing slab until homogeneous.

Apply the mixture thinly on the corresponding margins with a brush, i. e. on the fresh, dry cut and ground areas. The working time of **Ufi Gel SC** glazing is 1 min: 30 s. The setting time at room temperature is 10 min. Cooler conditions prolong the setting time. The setting time can be shortened to 3 min by using a pressure unit (40 °C - 45 °C).

Any earlier finishing should be carried out with great caution.

The **Ufi Gel SC** set contains 1 polishing disk REF 2049. Abrasive strips are available in varying grain sizes. The most suitable sizes are 80 (coarse) and 120 (medium). Use medium to low rotational speed (5,000 - 8,000 RPM) and apply only very light pressure when finishing the transition area. This prevents heat generation in the relining as well as excessive abrasion of the denture material.

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