


Universal Prosthetic Kit

Universal Prosthetic Kit

NeoBiotech Global HQ: Global Marketing

Neo
Biotech
Satisfaction to Dentists




Contents

I . Introduction

II . Features: General

III . Features: Components

IV . Summary



Introduction

Various types of driver in Universal Prosthetic Kit covers major types of prosthetic installation:
Straumann, Nobel Btiodcare, and NeoBiotech.



Features

✓ Simple


Only an Universal Prosthetic KIT is required for a hole procedure.

✓ Compatible

1.2 hex, 1.25 hex, SCS and Unigrip driver are compatible with othe brands



Features: General



Features: General

Neo Universal Prosthetic Kit Ver.0

Hex 1,2



15mm



10mm

Hex 1,25



15mm



10mm

Hex
Type
Driver

SCS



15mm



10mm

Unigrip




15mm



10mm

Other Connection
Type
Driver

Spare



Torque Ratchet Wrench

Assistant Tool

Driver holder

Neo
Biotech

Features: Components

Driver : Hex Type Connection

1.2 Hex


1.25 Hex

10mm

15mm

10mm

15mm



NeoBiotech Global HQ: Global Marketing




Neo Biotech
Satisfaction to Dentists

Features: Components

Driver : Other Type Connection

SGS Driver

Unigrip Driver



NeoBiotech Global HQ: Global Marketing

Neo Biotech
Satisfaction to Dentists

Features: Components

Compatibility List


1.2 Hex	1.25 Hex	SCS	Unigrip
			
<ul style="list-style-type: none"> ▪ Biomet 3i ▪ Dentsply-Friaden ▪ DIO ▪ Megazen ▪ Neobiotech ▪ Osstem ▪ Xive 	<ul style="list-style-type: none"> ▪ Astra ▪ Biohorizons ▪ Camlog ▪ Dentis ▪ Dentium ▪ Imlant Direct Legacy ▪ MIS ▪ Zimmer 	<ul style="list-style-type: none"> ▪ Straumann ▪ Implant Direct Replant 	<ul style="list-style-type: none"> ▪ Implant Direct Swishplait ▪ Neoss ▪ Nobel Biocare ▪ Southern Implants

Features: Components

Torque Ratchet

Installed on hex driver, Torque Ratchet tightens or looses prosthetic parts, and measures torque.


Direction tip:
pull & turn the arrow to change
direction of force applied



Features: Components

Driver Holder

Driver holder gives support drivers; hex, other than hex, or spares to intraorally tighten an abutment or abutment screw.



Thanks



Satisfaction to Dentists

NeoBiotech Global HQ: Global Marketing

