



TePe® Interdental Brush

Handle	Polypropylene (ISCC certified).
Tip	TePe G2™ neck (pink, orange, red and blue). Thermoplastic elastomer (TPE).
Cap	Biobased high-density polyethylene (HDPE), UV dry offset ink.
Filaments	Polyamide (PA).
Wire	Plastic coated stainless steel.

Packaging

1-pack	<i>Consumer pack, 1 pc/bag</i> <i>Retail carton, 150 pcs/carton</i> The bag is made of Polypropylene (PP).
2-pack	<i>Consumer pack, 2 pcs/bag</i> <i>Retail carton, 100 pcs/carton</i> The bag is made of Polypropylene (PP).
6-pack	<i>Consumer pack, 6 pcs/blister pack</i> <i>Retail carton, 10 packs/carton</i> The blister is made of polyethylene terephthalate (PET), partly recycled. The backing card is made of FSC certified cardboard. The tray in the retail box is made of polyethylene terephthalate (PET), partly recycled.
8-pack	<i>Consumer pack, 8 pcs/bag</i> <i>Retail carton, 10 bags/carton</i> The bag is made of biobased polyethylene.
25-pack	<i>Consumer pack, 25 brushes with cap*/bag</i> <i>Retail carton, 16 bags/carton</i> The bag is made of biobased polyethylene.

Retail and transport cartons for all packaging above:

Chlorine-free cardboard respectively corrugated cardboard. All from 100% recycled material, except for the retail carton for 6-pack which is from partly recycled material.

*Grey and black without cap.

tepe.com

TePe®