























STEP 1 CAVITY ACCESS & BURS

| LEXICON | CODE | ISO SIZE | mm | SYMBOL |
|--|--------------------|--|---|--|
| DIAMENDO  | A 0165 | |  |   534 C |
| ROOT CANAL DRILL  | A 0007 | 6x: N°1 Ø1.20mm, N°2 Ø1.40mm, N°3 Ø1.65mm, N°4 Ø1.90mm Assort. 2x N°1-N°2, 1x N°3-N°4 | |   800-1200 min ⁻¹ |
| GATES-GLIDDEN DRILL  | A 0008 | 6x: N°1 Ø 0.50mm, N°2 Ø 0.70mm, N°3 Ø 0.90mm, N°4 Ø1.10mm, N°5 Ø1.30mm, N°6 Ø1.50mm Assort. 1x: N°1, N°2, N°3, N°4, N°5, N°6 6x All sizes available in Sterile Packaging for lengths 28mm and 32mm. |  |   800-1200 min ⁻¹ |
| X-GATES LARGO® PEESO REAMER  | A 0008-7 A 0009 | 6x 6x: N°1 Ø 0.70mm, N°2 Ø 0.90mm, N°3 Ø1.10mm, N°4 Ø1.30mm, N°5 Ø1.50mm, N°6 Ø1.70mm Assort. 6x: N°1, N°2, N°3, N°4, N°5, N°6 All sizes available in Sterile Packaging. |  |   800-1200 min ⁻¹ |
| ENDO BURS | CODE | ISO SIZE | mm | SYMBOL |
| ENDO ACCESS BUR DR HOWARD MARTIN  | A 0164 | 1x: N°1 Ø1.1mm, N°2 Ø1.5mm, N°3 Ø1.8mm Assort. 1,2,3 |  |   534 C |
| ENDO-Z DR ZEKRYA  | E 0152 | 5x: 21mm, 25mm, 24.5mm |  |  |