

## Long Neck

E 0205

### Characteristics:

Long neck for enhanced visibility.

### Applications:

- Excavating carious dentine
- Widening and finishing cavity access
- Location of the root orifices



## Endo Access Bur

A 0164

### Characteristics:

Diamond plated conical bur with round tip (grit size 130µm).

### Applications:

- Opening and preparing the pulp cavity. Cavity access performed within one operation using a single bur on monoradicular teeth and shaping, on pluriradicular teeth, of a conical access adapted to the insertion of the Endo-Z and Diamendo burs



## Diamendo

A 0165

### Characteristics:

Diamond plated conical bur with non-cutting extremity to avoid perforating the floor of the pulp chamber or the root walls.

### Applications:

- Opening and preparing the pulp cavity (access to the root orifices)
- Rapid elimination of any interference before the insertion of canal instruments

