## AquaCem GI Luting Cement

## DeTrey® Zinc Zinc Phosphate Cement



\_Self-adhesive glass-ionomer luting cement that adheres to dental substances and restorations including:

- · Porcelain-fused to metal crowns or bridges
- Metal crowns or bridges, inlays and onlays
- Molar bands and post-core units
- Maryland bridges

The AquaCem Standard Package includes:

- 1 x Powder 30g
- 1 x Powder Measure
- 1 x Water Dispenser

## Features and Benefits

AquaCem is a light-yellow, translucent glass-ionomer luting material consisting of a blend of alumino-silicate glass and polyacrylic acid.

The powder is mixed with distilled water to produce a luting material which adheres to dentin and enamel producing tightly sealed cementations.