

























PRODUCT	Sycamore Interdental Wedges																																																													
REFERENCE	<table border="1"> <thead> <tr> <th colspan="2">Sycamore Interdental Wedges – Assortments</th> <th>REF/Part No.</th> </tr> </thead> <tbody> <tr> <td colspan="2">Wedge System Assortment + 1 Wedge Positioner Instrument – Pack of 1,000</td> <td>830</td> </tr> <tr> <td colspan="2">Sycamore Interdental Wedges – Assorted Set – Pack of 100</td> <td>823</td> </tr> <tr> <td colspan="2">Sycamore Interdental Wedges – Assorted Set – Pack of 400</td> <td>827</td> </tr> <tr> <td colspan="2">Sycamore Interdental Wedges – Wedges Dispenser – Assorted Set – Pack of 500</td> <td>826</td> </tr> <tr> <td colspan="2">Wedge Positioner Instrument</td> <td>820</td> </tr> </tbody> </table> <table border="1"> <thead> <tr> <th colspan="2">Sycamore Interdental Wedges – Refills</th> <th>REF/Part No.</th> </tr> </thead> <tbody> <tr> <td rowspan="2"></td> <td>Turquoise (Mini)</td> <td>Pack of 1,000 829/05</td> </tr> <tr> <td></td> <td>Pack of 100 822/05</td> </tr> <tr> <td rowspan="2"></td> <td>Orange</td> <td>Pack of 1,000 829/10</td> </tr> <tr> <td></td> <td>Pack of 100 822/10</td> </tr> <tr> <td rowspan="2"></td> <td>White</td> <td>Pack of 1,000 829/20</td> </tr> <tr> <td></td> <td>Pack of 100 822/20</td> </tr> <tr> <td rowspan="2"></td> <td>Green</td> <td>Pack of 1,000 829/30</td> </tr> <tr> <td></td> <td>Pack of 100 822/30</td> </tr> <tr> <td rowspan="2"></td> <td>Yellow</td> <td>Pack of 1,000 829/40</td> </tr> <tr> <td></td> <td>Pack of 100 822/40</td> </tr> <tr> <td rowspan="2"></td> <td>Blue</td> <td>Pack of 1,000 829/50</td> </tr> <tr> <td></td> <td>Pack of 100 822/50</td> </tr> <tr> <td rowspan="2"></td> <td>Red</td> <td>Pack of 1,000 829/60</td> </tr> <tr> <td></td> <td>Pack of 100 822/60</td> </tr> <tr> <td rowspan="2"></td> <td>White</td> <td>Pack of 1,000 829/70</td> </tr> <tr> <td></td> <td>Pack of 100 822/70</td> </tr> </tbody> </table>	Sycamore Interdental Wedges – Assortments		REF/Part No.	Wedge System Assortment + 1 Wedge Positioner Instrument – Pack of 1,000		830	Sycamore Interdental Wedges – Assorted Set – Pack of 100		823	Sycamore Interdental Wedges – Assorted Set – Pack of 400		827	Sycamore Interdental Wedges – Wedges Dispenser – Assorted Set – Pack of 500		826	Wedge Positioner Instrument		820	Sycamore Interdental Wedges – Refills		REF/Part No.		Turquoise (Mini)	Pack of 1,000 829/05		Pack of 100 822/05		Orange	Pack of 1,000 829/10		Pack of 100 822/10		White	Pack of 1,000 829/20		Pack of 100 822/20		Green	Pack of 1,000 829/30		Pack of 100 822/30		Yellow	Pack of 1,000 829/40		Pack of 100 822/40		Blue	Pack of 1,000 829/50		Pack of 100 822/50		Red	Pack of 1,000 829/60		Pack of 100 822/60		White	Pack of 1,000 829/70		Pack of 100 822/70
Sycamore Interdental Wedges – Assortments		REF/Part No.																																																												
Wedge System Assortment + 1 Wedge Positioner Instrument – Pack of 1,000		830																																																												
Sycamore Interdental Wedges – Assorted Set – Pack of 100		823																																																												
Sycamore Interdental Wedges – Assorted Set – Pack of 400		827																																																												
Sycamore Interdental Wedges – Wedges Dispenser – Assorted Set – Pack of 500		826																																																												
Wedge Positioner Instrument		820																																																												
Sycamore Interdental Wedges – Refills		REF/Part No.																																																												
	Turquoise (Mini)	Pack of 1,000 829/05																																																												
		Pack of 100 822/05																																																												
	Orange	Pack of 1,000 829/10																																																												
		Pack of 100 822/10																																																												
	White	Pack of 1,000 829/20																																																												
		Pack of 100 822/20																																																												
	Green	Pack of 1,000 829/30																																																												
		Pack of 100 822/30																																																												
	Yellow	Pack of 1,000 829/40																																																												
		Pack of 100 822/40																																																												
	Blue	Pack of 1,000 829/50																																																												
		Pack of 100 822/50																																																												
	Red	Pack of 1,000 829/60																																																												
		Pack of 100 822/60																																																												
	White	Pack of 1,000 829/70																																																												
		Pack of 100 822/70																																																												
COMPANY/ Manufacturer	KerrHawe SA																																																													
CLASSIFICATION	CLASS I MEDICAL DEVICE																																																													
PRODUCT DESCRIPTION	Interdental Wedges with anatomical shape. The concave sides adapt perfectly to the interdental contour. The square end allows easy insertion of the wedge and the upturned tip prevents damage to the papillae.																																																													
SPECIFICATIONS	<p>The wedges are made of splinter-free sycamore wood (technical specifications in the table below)</p> <p>Materials for item 820: stainless steel AISI -304 and AISI 420</p>																																																													

	Articles	Color (color material)	Length	Base Section
	829/05 822/05	Turquoise (E-131)	12 mm	1.2 mm x 1.6 mm
	829/10 822/10	Orange (E-110)	13 mm	1.4 mm x 2 mm
	829/20 822/20	White (No color)	13 mm	2 mm x 2 mm
	829/30 822/30	Green (E-131 and E-104)	13 mm	2 mm x 2 mm
	829/40 822/40	Yellow (E-104)	13 mm	2 mm x 2 mm
	829/50 822/50	Blue (E-131)	15 mm	2 mm x 2 mm
	829/60 822/60	Red (E-122)	16 mm	2.5 mm x 2.5 mm
	829/70 822/70	White (No color)	18 mm	2.5 mm x 2.5 mm
PACKAGING	Transparent packaging 5.8 cm width 8.6 cm length, 2.9 cm height			
SHELF LIFE	1800 days			
INTENDED OF USE	<ul style="list-style-type: none"> • Adaptation of the matrix along the entire cervical margin, • Separation during the filling phase (contact points) • Stability of matrix by means of continuous pressure zone 			
SECURITY	Only profesional use			
TRANSPORT AND STORAGE CONDITIONS	Ambient temperature			