


<b>PRODUCT</b>	Hawe Tofflemire Contoured Matrices										
<b>REFERENCE</b>	<table border="1"> <thead> <tr> <th colspan="3">Hawe Tofflemire Contoured Matrices</th> </tr> </thead> <tbody> <tr> <td>Refill pack of 30</td> <td>Hawe Tofflemire Contoured Matrices, 0.050mm thin</td> <td>1001 C</td> </tr> <tr> <td></td> <td>Hawe Tofflemire Contoured Matrices, 0.038mm thin</td> <td>1101 C</td> </tr> </tbody> </table>			Hawe Tofflemire Contoured Matrices			Refill pack of 30	Hawe Tofflemire Contoured Matrices, 0.050mm thin	1001 C		Hawe Tofflemire Contoured Matrices, 0.038mm thin
Hawe Tofflemire Contoured Matrices											
Refill pack of 30	Hawe Tofflemire Contoured Matrices, 0.050mm thin	1001 C									
	Hawe Tofflemire Contoured Matrices, 0.038mm thin	1101 C									
<b>COMPANY/ Manufacturer</b>	KerrHawe SA										
<b>CLASSIFICATION</b>	CLASS I MEDICAL DEVICE										
<b>PRODUCT DESCRIPTION</b>	Countoured matrices										
<b>SPECIFICATIONS</b>	<p>1001 C:</p> <ul style="list-style-type: none"> <li>• Material: Stainless steel 1.4301</li> <li>• Overall dimensions: 63.4 mm x 13.01 mm</li> <li>• Thickness: 0.050 mm</li> </ul> <p>1101 C</p> <ul style="list-style-type: none"> <li>• Material: Stainless steel 1.4310</li> <li>• Overall dimensions: 63.4 mm x 13.01 mm</li> <li>• Thickness: 0.038 mm</li> </ul>										
<b>PACKAGING</b>	Refill pack of 30 pcs: 5.8 cm width 8.6 cm length, 2.9 cm height										
<b>SHELF LIFE</b>	1800 days										
<b>INTENDED OF USE</b>	Use for both molar and premolar restorations										
<b>SECURITY</b>	Only professional use										
<b>TRANSPORT AND STORAGE CONDITIONS</b>	Ambient temperature										